English	3
Bulgarian (Български език)	7
Chinese (Simplified) (简体中文).....	11
Croatian (Hrvatski).....	15
Czech (Čeština).....	19
Danish (Dansk).....	23
Dutch (Nederlands)	27
Estonian (Eesti)	31
Finnish (Suomi)	35
French (Français)	39
German (Deutsch).....	43
Greek (Ελληνικά).....	47
Hungarian (Magyar).....	51
Indonesian (Bahasa Indonesia).....	55
Italian (Italiano)	59
Japanese (日本語).....	63
Korean (한국어).....	67
Latvian (Latviešu)	71
Lithuanian (Lietuvių)	75
Norwegian (Norsk).....	79
Polish (Polski).....	83
Portuguese (Portugal) (Português (Portugal)).....	87
Romanian (Română)	91
Russian (Русский)	95
Serbian (Srpski).....	99
Slovak (Slovensky)	103
Spanish (Español)	107
Swedish (Svenska).....	111
Turkish (Türkçe)	115
Vietnamese (Tiếng Việt)	119

Blank Page

GE Healthcare

PLEASE READ THIS FIRST

Advantage Workstation*

DICOM VIEWER TROUBLESHOOTING GUIDE

5447852-199

Revision 1

© 2011 General Electric Company.

Legal Notices

*GE, the GE Monogram, GE Healthcare, Centricity, and Advantage Workstation are trademarks of General Electric Company.

Windows, Windows Vista and the Windows logo are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

All other products and company names mentioned herein may be the trademarks of their respective owners.

Advantage Workstation – DICOM Viewer Troubleshooting Guide is an independent publication and is not affiliated with, nor has it been authorized, sponsored, or otherwise approved by Microsoft Corporation.

www.gehealthcare.com

DICOM Viewer Troubleshooting Guide

This guide is applicable only if DICOM Viewer is used on supported Windows® operating systems.

Note: Centricity* DICOM Viewer 3.1.4 is supported on the following Windows® operating systems:

- 32-bit Windows® XP (including Service Pack 2 and Service Pack 3)
- 32-bit Windows Vista® (including Service Pack 1 and Service Pack 2)
- 32-bit Windows® 7

Problem	Solution
DICOM Viewer does not run automatically.	<p>Disable the 'Turn off Autoplay' feature on the Windows® operating system as follows:</p> <p><i>On Windows® XP:</i></p> <ol style="list-style-type: none"> 1. Press <Windows key + R> to open Run, type <code>gpedit.msc</code>, and then click (OK). 2. In the Group Policy window, under Computer Configuration, expand Administrative Templates, and click System. 3. In the System/Setting pane, double-click Turn off Autoplay. 4. Select the Disabled check box, click (Apply), and then click (OK). <p><i>On Windows Vista®/Windows® 7:</i></p> <ol style="list-style-type: none"> 1. Press <Windows key + R> to open Run, type <code>gpedit.msc</code>, and then click (OK). 2. In the Local Group Policy Editor window, under Computer Configuration, expand Administrative Templates, expand Windows Components, and then click AutoPlay Policies. 3. In the Setting pane, double-click Turn off Autoplay. 4. Select the Disabled check box, click (Apply), and then click (OK). <p>Note: Gpedit.msc may not be available in Windows® XP, Windows Vista® and Windows® 7 Home editions.</p>
DICOM Viewer gives an error about compression.	<p>You must recreate data and set the CD/DVD profile to General Purpose on the AW as follows:</p> <ol style="list-style-type: none"> 1. Open Media Composer. 2. Click the (Preferences) button. 3. Set profile to General Purpose.
DICOM Viewer installation fails.	<p>On the station with Windows® operating system, delete the temporary 'dicomviewer' directory located under: D:\Documents and Settings\<UserID>\Local Settings\Temp\</p> <p>Note that this directory structure may also be located on the C: drive.</p>

If the security settings prevent performing the specified steps or you encounter any other problems, please contact the site's Laptop/PC or IT admin support for assistance.

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